

Stride One in Gait & Mobility (Re)training:

Improving Recovery After TKA, THA, ACL & Achilles Tendon Surgery



Following joint surgery, the majority of patients report decreased pain and successful long-term results, but **recovery is unpredictable**, and most patients continue to exhibit muscle weakness in their lower limbs and functional limitations in comparison to similarly aged control individuals¹. One month after a total knee replacement for example, muscle function was shown to be lowered by 20%-25% but one year later, it is **still lower than in healthy adults**, with reports of 18% slower walking speed and 51% slower stair climbing speed^{2,3}.

Habitual **movement compensations**, such as decreased peak surgical knee extension moments persist years after TKA, are linked to poorer recovery, and may influence contralateral osteoarthritis progression⁴. Gait compensation (e.g. **asymmetry**) remains similar from 3- to 6-months during a task requiring greater knee demand compared to overground walking post-TKA, despite improvements in self-report surveys⁵.

Surgeons experience **inconsistencies in patient satisfaction** prediction. Patients often do not reach functional levels of healthy controls and present with persistent gait deviations^{6,7,8}.

Evidence suggests that **biofeedback in early postoperative rehabilitation** after TKA is effective in improving gait symmetry, reducing pain and increasing activity level⁹. Weight-bearing biofeedback training following TKA results in increased knee extension moments during gait and improved the times measured for the sit-to-stand tests¹. The most common types of feedback used are visual feedback followed by audio feedback, with one study mentioning that audio was preferred by the patients as it was easier to handle. Overall, the large majority of feedback methods demonstrated a potential value for improving mobility and decreasing pain⁶.

References:

1. Dutta S, Ambade R, Wankhade D, Agrawal P. Rehabilitation Techniques Before and After Total Knee Arthroplasty for a Better Quality of Life. *Cureus*. 2024 Feb 25;16(2):e54877
2. Petterson SC, Mizner RL, Stevens JE, Raisis L, Bodenstab A, Newcomb W, Snyder-Mackler L. Improved function from progressive strengthening interventions after total knee arthroplasty: a randomized clinical trial with an imbedded prospective cohort. *Arthritis Rheum*. 2009 61:174–183.
3. Bade MJ, Kohrt WM, Stevens-Lapsley JE. Outcomes before and after total knee arthroplasty compared to healthy adults. *J Orthop Sports Phys Ther*. 2010 40:559–567.
4. Christiansen CL, Bade MJ, Davidson BS, Dayton MR, Stevens-Lapsley JE. Effects of Weight-Bearing Biofeedback Training on Functional Movement Patterns Following Total Knee Arthroplasty: A Randomized Controlled Trial. *J Orthop Sports Phys Ther*. 2015 Sep;45(9):647-55.
5. Christensen JC, Pelt CE, Bo Foreman K, LaStayo PC, Anderson AE, Gililand JM, Mizner RL. Longitudinal study of knee load avoidant movement behavior after total knee arthroplasty with recommendations for future retraining interventions. *Knee*. 2021 Jun;30:90-99. doi: 10.1016/j.knee.2021.03.014.
6. Howell, D. (2025). Total Knee Joint Replacement – Are You Satisfied? Gait Analysis & Training Can Provide Post-Op Improvements.
7. Al-Amri, M. N., et al. (2020). Patient-Reported Outcomes Following Total Knee Replacement in Patients <65 Years of Age—A Systematic Review and Meta-Analysis. *Journal of Clinical Medicine*, 9(10), 3150.
8. Al-Amri, M. N., et al. (2024). Inconsistencies in survey-based (eg, questionnaire) data collection across biomedical, clinical, behavioral, and social sciences pose challenges to research reproducibility
9. Pfeufer D, Gililand J, Böcker W, Kammerlander C, Anderson M, Krähenbühl N, Pelt C. Training with biofeedback devices improves clinical outcome compared to usual care in patients with unilateral TKA: a systematic review. *Knee Surg Sports Traumatol Arthrosc*. 2019 May;27(5):1611-1620.

Improving Recovery After TKA, THA, ACL & Achilles Tendon Surgery



In people with TKA sometimes partial weight bearing is prescribed per post-operative rehabilitation protocol. However, patients are **unable to perform a correct gait with low weight bearing instructions** without a smart insole with biofeedback and measuring system¹⁰. On the other hand, using biofeedback from a wireless insole sensor significantly improves the maximum peak load in patients with prescribed unrestricted weight bearing¹¹. This suggests improved **compliance and confidence with weight bearing instructions**, a precondition for gait and mobility training.

Currently, most postoperative rehabilitation from TKA occurs in the outpatient setting or at home. This leads clinicians to rely on reports from outpatient physical therapy or subjective methods such as patient-reported outcomes measures (PROMs) to measure recovery. However, previous studies have outlined **concerns over the standardization of PROMs as the only measure of recovery**. For instance, when determining if patients accurately report distance walked, the mean magnitude of reporting error is 69% preoperatively and 93% postoperatively.

Providers should exercise caution when interpreting patient reported activity levels before and after TKA/THA¹². Moreover, clinic visits to evaluate recovery have often been reduced to 4–6 weeks postoperatively, increasing the possibility for important landmarks of recovery to be missed. This raises the potential value of objective activity measurement and the consistent contact between patients and surgical teams that wearable devices offer following TKA^{13,14}.

References:

10. Castellarin, Gianluca & Merlini, Michele & Bettinelli, Giulia & Riso, Raffaella & Bori, Edoardo & Innocenti, Bernardo. (2022). Effect of an Innovative Biofeedback Insole on Patient Rehabilitation after Total Knee Arthroplasty. Applied Sciences. 12. 2456. 10.3390/app12052456.
11. M. Raaben, S. Redzwan, R. Augustine, T.J. Blokhuis. Real-time biofeedback improves rehabilitation after hip fractures in elderly: an international, multicenter randomised controlled trial (COMFORT).
12. Vaughn NH, Dunkleberger MF, Mason MW. Individual patient-reported activity levels before and after joint arthroplasty are neither accurate nor reproducible. Clinical Orthopaedics and Related Research 2019477536–544.
13. Luna IE, Kehlet H, Peterson B, Wede HR, Hoevsgaard SJ, Aasvang EK. Early patient-reported outcomes versus objective function after total hip and knee arthroplasty: a prospective cohort study. Bone and Joint Journal 201799-B1167–1175.
14. Bandholm T, Wainwright TW, Kehlet H. Rehabilitation strategies for optimisation of functional recovery after major joint replacement. Journal of Experimental Orthopaedics 20185 44.

Improving Recovery After TKA, THA, ACL & Achilles Tendon Surgery



Assessment of gait quality using spatial-temporal gait parameters (STGP) has proven applications in measuring fall risk, monitoring postoperative changes in mobility, and optimizing treatment plans. Additionally, knee biomechanics early after knee replacement surgery predict abnormal gait patterns 12 months postoperatively^{15,16,17}.

Patients receiving feedback on gait parameters such as step counts have significantly higher mean daily step count by 43% in week 1, 33% in week 2, 21% in week 6, and 17% at 6 months. Commercial activity trackers with step count feedback are associated with **higher activity levels** after TKA/THA and may be a useful adjunct after surgery¹⁸.

Post-operative care costs of lower limb arthroplasty account for 36% of the total episode of care costs. A platform for clinicians that allows them to view patient engagement, activity levels, and gait parameters demonstrated non-inferiority to traditional care models (no significant difference in 90-day mean flexion, 90-day mean single leg stance time, 90-day mean timed up and go time) and potential to **decrease postoperative costs** while improving patient engagement and communication with providers^{19,20}. Mehta et al. showed that goal setting and connection to care teams significantly reduced rehospitalizations²¹.

References:

15. Hausdorff JM. Gait variability: methods, modeling and meaning example of increased stride time variability in elderly fallers quantification of stride-to-stride fluctuations. *Journal of NeuroEngineering and Rehabilitation* 2005;9:1–9.
16. Weiss A, Brozgov M, Dorfman M, Herman T, Shema S, Giladi N, Hausdorff JM. Does the evaluation of gait quality during daily life provide insight into fall risk? A novel approach using 3-day accelerometer recordings. *Neurorehabilitation and Neural Repair* 2013;27:742–752.
17. Levinger P, Menz HB, Morrow AD, Perrott MA, Bartlett JR, Feller JA, Bergman NB. Knee biomechanics early after knee replacement surgery predict abnormal gait patterns 12 months postoperatively. *Journal of Orthopaedic Research* 2012;30:371–376.
18. Van der Walt N, Salmon LJ, Gooden B, Lyons MC, O'Sullivan M, Martina K, Pinczewski LA, Roe JP. Feedback from activity trackers improves daily step count after knee and hip arthroplasty: a randomized controlled trial. *Journal of Arthroplasty* 2018;33:422–428.
19. Bozic KJ, Ward L, Vail TP, Maze M. Bundled payments in total joint arthroplasty: targeting opportunities for quality improvement and cost reduction. *Clinical Orthopaedics and Related Research* 2014;472:188–193.
20. Crawford DA, Duwelius PJ, Sneller MA, Morris MJ, Hurst JM, Berend KR, Lombardi AV. 2021 Mark Coventry Award: use of a smartphone-based care platform after primary partial and total knee arthroplasty: a prospective randomized controlled trial. *Bone and Joint Journal* 2021;103-B.
21. Mehta SJ, Hume E, Troxel AB, Reitz C, Norton L, Lacko H, McDonald C, Freeman J, Marcus N, Volpp KG et al. Effect of remote monitoring on discharge to home, return to activity, and rehospitalization after hip and knee arthroplasty: a randomized clinical trial. *JAMA Network Open* 2020;3:e2028328.

Improving Recovery After TKA, THA, ACL & Achilles Tendon Surgery



What's next?

Reliable digital and real-time data on patients' gait quality that are accessible to both the surgeon and the patient may help with **steering the recovery process and improve patient satisfaction**^{22,23,24,25,26,27,28}.

Real-time gait analysis and biofeedback to immediately identify and correct asymmetries, weight-bearing issues, and compensatory patterns may result in data-driven adjustments to exercises, optimizing ROM and strength gains for faster, more effective recovery. Visual feedback and progress tracking increase motivation and adherence, leading to better long-term outcomes and satisfaction. Objective gait data for precise monitoring of patient progress allow for timely and evidence-based adjustments to rehabilitation protocols.

Wearables providing real-time gait feedback have the potential to **supplement or completely replace conventional rehabilitation regimens** by effectively tracking physical activity and enhancing patient participation after total knee arthroplasty. Intervention components, including step objectives, app-based patient interaction platforms, and patient-specific rehabilitation benchmarks, may help programs work better²⁹. Self-directed rehabilitation with a wearable device linked to a smartphone app showed non-inferiority to standard rehabilitation following TKA^{30,31}. A significant better motor performance (gait symmetry and single stance support) has been found when applying visual feedback gait training after TKA compared with conventional physical therapy³².

References:

22. Kayani B, Konan S, et al. Five-year outcomes of minimally invasive total knee replacement. Bone Joint J 2019.
23. Wernecke GC, et al. Prevalence of knee stiffness after TKR and strategies to minimize it: J Arthroplasty 2015.
24. Pizzari T et al. Quadriceps activation after knee surgery: J Orthop Phys Ther.
25. Maniar RN, et al. Early rehabilitation after TKR: protocols and gait.
26. Smith TO, et al. Return to driving after hip and knee arthroplasty.
27. Wickstrom E, et al. Adherence and motivation during post-op gait rehabilitation.
28. McClelland JA, et al. Patient adherence and function post TKR.
29. Constantinescu D, Pavlis W, Rizzo M, Vanden Berge D, Barnhill S, Hernandez VH. The role of commercially available smartphone apps and wearable devices in monitoring patients after total knee arthroplasty: a systematic review. EFORT Open Rev. 2022 Jul 5;7(7):481-490.
30. Crawford DA, Duwelius PJ, Sneller MA, Morris MJ, Hurst JM, Berend KR, Lombardi AV. 2021 Mark Coventry Award: use of a smartphone-based care platform after primary partial and total knee arthroplasty: a prospective randomized controlled trial. Bone and Joint Journal 2021103-B (6 Supplement) 3–12.
31. Tripuraneni KR, Foran JRH, Munson NR, Racca NE, Carothers JT. A smartwatch paired with a mobile application provides postoperative self-directed rehabilitation without compromising total knee arthroplasty outcomes: a randomized controlled trial. Journal of Arthroplasty 2021363888–3893.
32. Carozzo S, Vatrano M, Coschignano F, Battaglia R, Calabrò RS, Pignolo L, Contrada M, Tonin P, Cerasa A, Demeco A. Efficacy of Visual Feedback Training for Motor Recovery in Post-Operative Subjects with Knee Replacement: A Randomized Controlled Trial. J Clin Med. 2022 Dec 11;11(24):7355.

Improving Recovery After TKA, THA, ACL & Achilles Tendon Surgery



Wearable devices and smartphone apps are able to **increase physical activity following TKA**. Van der Walt et al. showed setting daily step goals and using wearable devices that provide step count feedback can significantly increase daily step count up to 6 months postoperatively¹⁷.

Preoperative step count, change in step count, and postoperative step count can be used to **predict patient recovery**.

These findings suggest value in leveraging passively collected step data from patients before surgery. While more research needs to be done to set data-driven, individualized benchmarks, wearable data can be a clinically viable alternative to predict recovery outcomes and identify complications earlier^{33,34}

Join our pilot program

Scan* the QR code to learn more and apply.



*Open your smartphone camera and point it at the QR code to get started.

Ceriter Nederland BV,
François de Veyestraat 8b
6211AB Maastricht
The Netherlands

P: 32 (0)89 39 21 64
E: info@ceriter.com

Disclaimer:

Stride One is a Class I medical device certified under the Medical Device Directive 93/42/EEC. It conforms to applicable EU safety and performance requirements and bears the CE mark. This certification applies solely to the product's current intended use and classification under MDD. Stride One is intended to support post-operative rehabilitation following lower-limb orthopedic procedures. It is not intended to replace clinical judgment, professional physiotherapy, or prescribed rehabilitation protocols. Use should be advised by a healthcare professional. Clinical decisions should be based on the full clinical context, including patient-specific factors and professional expertise. Patient selection, rehabilitation protocols, and monitoring remain the sole responsibility of the treating healthcare provider. The device is currently undergoing further clinical validation. Ceriter makes no claims regarding treatment outcomes or specific levels of efficacy beyond the scope of its certified use.

References:

33. Bini SA, Shah RF, Bendich I, Patterson JT, Hwang KM, Zaid MB. Machine learning algorithms can use wearable sensor data to accurately predict six-week patient-reported outcome scores following joint replacement in a prospective trial. *Journal of Arthroplasty* 2019;34:2242–2247.
34. Patterson JT, Wu HH, Chung CC, Bendich I, Barry JJ, Bini SA. Wearable activity sensors and early pain after total joint arthroplasty. *Arthroplasty Today* 2020;668–70.